



# yashima



## Brinell Type Shop Microscope



A simple highly accurate instrument to permit exact measurements. Supplied with adjustable penlight, using standard batteries. For testing surface of Metal, Fibre, Precisely Worked Parts and so on. This shop microscope is well built to perform under all conditions.

Type	Magnification	Scale
A	20 ×	0.05 $\frac{m}{m}$
B	60 ×	0.02 $\frac{m}{m}$
C	100 ×	0.01 $\frac{m}{m}$

Special order	Magnification	Scale
	10 ×	0.10 $\frac{m}{m}$
	30 ×	0.05 $\frac{m}{m}$
	40 ×	0.025 $\frac{m}{m}$
	50 ×	0.02 $\frac{m}{m}$

"inch" scale also acceptable by special order.

**YASHIMA OPPICEL CO.,**